

Microscope Checklist

Did you.....

- 1. turn off the lamp?
- 2. put objective lens at its lowest setting (usually 4x)?
- 3. clean & return slide & cover slip?
- 4. set diaphragm at 3?
- 5. tuck electrical cord away neatly?
- 6. clean any dirty lens and the stage?
- 7. cover the microscope?

NEVER

1. Place microscope near the edge of the table.
2. Use coarse adjustment knob when using high power objective lens (40x and above).
3. Move stage toward the objective lens when focusing.